

Datasheet

Hsp90b1 polyclonal antibody (RPE)

gp96

Catalog Number: PAB28962**Regulatory Status:** For research use only (RUO)**Product Description:** Rabbit polyclonal antibody raised against synthetic peptide of Hsp90b1.**Immunogen:** A synthetic peptide corresponding to amino acids at C-terminus of mouse Hsp90b1.**Host:** Rabbit**Reactivity:** Bovine, Human, Mouse, Rat**Applications:** ICC, IHC, IP, WB

(See our web site product page for detailed applications information)

Protocols: See our web site at<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols**Specificity:** Detects ~94 kDa.**Form:** Liquid**Conjugation:** RPE**Purification:** Peptide Affinity Purified**Recommend Usage:** Immunocytochemistry

Immunohistochemistry

Immunoprecipitation (1:80)

Western blot (1:1000)

The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH 7.4 (50% glycerol, 0.09% sodium azide)**Storage Instruction:** Store at 4°C.

Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 22027**Gene Symbol:** Hsp90b1**Gene Alias:** ERp99, GRP94, TA-3, Targ2, Tra-1, Tra1,