

Datasheet

KDEL polyclonal antibody (PerCP)

Catalog Number: PAB28979

Regulatory Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against synthetic peptide of KDEL.

Immunogen: A synthetic peptide corresponding to KDEL.

Host: Rabbit

Reactivity: Human, Mouse, Rat

Applications: IHC, WB

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: Recognizes KDEL proteins, Hsp90b1, Hspa5, PDI and Calr. It may also see PDIA3 and Erp29.

Form: Liquid

Conjugation: PerCP

Purification: Affinity purification

Recommend Usage: Immunohistochemistry

Western blot (1:1000)

The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH 7.2 (50% glycerol, 0.09% sodium azide)

Storage Instruction: Store at 4°C.

Aliquot to avoid repeated freezing and thawing.