

Datasheet

Sod2 polyclonal antibody (PerCP)

Entrez GeneID: 24787**Catalog Number:** PAB28997**Gene Symbol:** Sod2**Regulatory Status:** For research use only (RUO)**Gene Alias:** -**Product Description:** Rabbit polyclonal antibody raised against Recombinant of Sod2.**Immunogen:** Recombinant protein corresponding rat Sod2.**Host:** Rabbit**Reactivity:** Bovine, Chicken, Dog, Fruit fly, Guinea pig, Hamster, Human, Monkey, Mouse, Pig, Rabbit, Rat, Sheep, Xenopus**Applications:** ELISA, IF, IHC, IP, WB

(See our web site product page for detailed applications information)

Protocols: See our web site at<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols**Specificity:** Detects a ~25 kDa protein corresponding to the molecular mass of Mn superoxide dismutase (SOD) on SDS PAGE immunoblots.**Form:** Liquid**Conjugation:** PerCP**Purification:** Antigen affinity purification**Recommend Usage:** ELISA

Immunofluorescence

Immunohistochemistry

Immunoprecipitation

Western blot (0.5 ug/ml)

The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH 7.4 (50% glycerol, 0.09% sodium azide)**Storage Instruction:** Store at 4°C.

Aliquot to avoid repeated freezing and thawing.