

Datasheet

RAB5A polyclonal antibody (FITC)

Catalog Number: PAB29012

Regulatory Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against synthetic peptide of RAB5A.

Immunogen: A synthetic peptide conjugated with KLH corresponding to human RAB5A.

Host: Rabbit

Reactivity: Bovine, Human, Monkey, Mouse, Rat

Applications: IF, WB

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: The antibody recognizes 26 kDa Rab5 of the above species origins.

Form: Liquid

Conjugation: FITC

Purification: Affinity purification

Recommend Usage: Immunofluorescence

Western blot (1:1000)

The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH 7.4 (50% glycerol, 0.09% sodium azide)

Storage Instruction: Store at 4°C.

Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 5868

Gene Symbol: RAB5A

Gene Alias: RAB5