

## Datasheet

### RAB5A polyclonal antibody (PerCP)

**Catalog Number:** PAB29015

**Regulatory Status:** For research use only (RUO)

**Product Description:** Rabbit polyclonal antibody raised against synthetic peptide of RAB5A.

**Immunogen:** A synthetic peptide conjugated with KLH corresponding to human RAB5A.

**Host:** Rabbit

**Reactivity:** Bovine, Human, Monkey, Mouse, Rat

**Applications:** IF, WB

(See our web site product page for detailed applications information)

**Protocols:** See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Specificity:** The antibody recognizes 26 kDa Rab5 of the above species origins.

**Form:** Liquid

**Conjugation:** PerCP

**Purification:** Affinity purification

**Recommend Usage:** Immunofluorescence

Western blot (1:1000)

The optimal working dilution should be determined by the end user.

**Storage Buffer:** In PBS, pH 7.4 (50% glycerol, 0.09% sodium azide)

**Storage Instruction:** Store at -20°C.

Aliquot to avoid repeated freezing and thawing.

**Entrez GeneID:** 5868

**Gene Symbol:** RAB5A

**Gene Alias:** RAB5