

Datasheet

Slc9a3 polyclonal antibody (PerCP)

Catalog Number: PAB29033

Regulatory Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against synthetic peptide of Slc9a3.

Immunogen: A synthetic peptide corresponding to amino acids 809-831 of rat Slc9a3.

Host: Rabbit

Applications: ICC, WB

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: Detects ~84 kDa.

Form: Liquid

Conjugation: PerCP

Purification: Affinity purification

Recommend Usage: Immunocytochemistry

Western blot (1:1000)

The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH 7.4 (50% glycerol, 0.09% sodium azide)

Storage Instruction: Store at 4°C.

Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 24784

Gene Symbol: Slc9a3

Gene Alias: Nhe3