

Datasheet

Actin polyclonal antibody

Catalog Number: PAB29054

Regulatory Status: For research use only (RUO)

Product Description: Chicken polyclonal antibody raised against purified Rabbit Actin.

Immunogen: Purified rabbit muscle actin.

Host: Chicken

Reactivity: Human, Mouse, Rabbit, Rat

Applications: ICC, WB

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Purification: Antigen affinity purification

Isotype: IgY

Recommend Usage: Immunocytochemistry

Western Blot(1:2500-1:5000)

The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH 7.2 (50% glycerol, 0.02% sodium azide)

Storage Instruction: Store at -20°C and avoid from light.

Aliquot to avoid repeated freezing and thawing.