

## Datasheet

### Beta-Galactosidase polyclonal antibody

**Catalog Number:** PAB29055

**Regulatory Status:** For research use only (RUO)

**Product Description:** Chicken polyclonal antibody raised against purified beta-galactosidase.

**Immunogen:** Purified beta-galactosidase protein (*lac Z* from *E. coli*)

**Host:** Chicken

**Applications:** ICC, IHC, WB  
(See our web site product page for detailed applications information)

**Protocols:** See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Form:** Liquid

**Purification:** Antigen affinity purification

**Isotype:** IgY

#### Recommend

**Usage:** Immunocytochemistry(1:2500-1:5000)

Immunohistochemistry(1:2500-1:5000)

Western Blot(1:5000-1:10000)

The optimal working dilution should be determined by the end user.

**Storage Buffer:** In PBS, pH 7.2 (50% glycerol, 0.01% thimerosal)

**Storage Instruction:** Store at 4°C. Store at -20°C for long term storage after dilution with glycerol and avoid from light.

Aliquot to avoid repeated freezing and thawing.