

Datasheet

FLAG polyclonal antibody

Catalog Number: PAB29056

Regulatory Status: For research use only (RUO)

Product Description: Chicken polyclonal antibody raised against recombinant FLAG peptide.

Immunogen: A KLH-conjugated synthetic peptide corresponding to Flag tag emulsified in Freund's adjuvant.

Sequence: DYKDDDDK

Host: Chicken

Applications: ELISA, ICC, IHC, IP, WB
(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Purification: Antigen affinity purification

Isotype: IgY

Recommend Usage: ELISA
Immunocytochemistry(1:500)
Immunohistochemistry
Immunoprecipitation(1:500)
Western Blot(1:1000)

The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH 7.2 (0.02% sodium azide)

Storage Instruction: Store at 4°C and avoid from light.
For long term storage, dilute the antibody with sterile glycerol and store it at -20°C.
Aliquot to avoid repeated freezing and thawing.