

Datasheet

HA polyclonal antibody (FITC)

Catalog Number: PAB29058

Regulatory Status: For research use only (RUO)

Product Description: Chicken polyclonal antibody raised against recombinant HA peptide.

Immunogen: A KLH-conjugated synthetic peptide corresponding to HA tag emulsified in Freund's adjuvant.

Sequence: YPYDVPDYA

Host: Chicken

Applications: ICC, IHC, WB
(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Conjugation: FITC

Purification: Antigen affinity purification

Isotype: IgY

Recommend Usage: Immunocytochemistry
Immunohistochemistry(1:250)
Western Blot
The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH 7.2 (0.2% BSA, 0.02% sodium azide)

Storage Instruction: Store at 4°C and avoid from light.
For long term storage, dilute the antibody with sterile glycerol and store it at -20°C.
Aliquot to avoid repeated freezing and thawing.