

Datasheet

HIS polyclonal antibody

Catalog Number: PAB29062

Regulatory Status: For research use only (RUO)

Product Description: Chicken polyclonal antibody raised against recombinant HIS peptide.

Immunogen: A KLH-conjugated synthetic peptide corresponding to 6X-Histidine tag emulsified in Freund's adjuvant.

Sequence: HHHHHH

Host: Chicken

Applications: ChIP, IP
(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Purification: Antigen affinity purification

Isotype: IgY

Recommend Usage: Chromatin Immunoprecipitation(1:5000-1:10000)
Immunoprecipitation(1:5000-1:10000)
The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH 7.2 (50% glycerol, 0.01% thimerosal)

Storage Instruction: Store at -20°C and avoid from light.
Aliquot to avoid repeated freezing and thawing.