

## Datasheet

### HIS polyclonal antibody

**Catalog Number:** PAB29063

**Regulatory Status:** For research use only (RUO)

**Product Description:** Chicken polyclonal antibody raised against recombinant HIS peptide.

**Immunogen:** A KLH-conjugated synthetic peptide corresponding to 6X-Histidine tag emulsified in Freund's adjuvant.

**Sequence:** HHHHHH

**Host:** Chicken

**Applications:** WB

(See our web site product page for detailed applications information)

**Protocols:** See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Form:** Liquid

**Purification:** Antigen affinity purification

**Isotype:** IgY

**Recommend Usage:** Western Blot(1:5000-1:10000)

The optimal working dilution should be determined by the end user.

**Storage Buffer:** In PBS, pH 7.2 (50% glycerol, 0.01% thimerosal)

**Storage Instruction:** Store at -20°C and avoid from light.

Aliquot to avoid repeated freezing and thawing.