

## Datasheet

### CHAT polyclonal antibody

**Catalog Number:** PAB29065

**Regulatory Status:** For research use only (RUO)

**Product Description:** Chicken polyclonal antibody raised against recombinant Human CHAT.

**Immunogen:** A KLH-conjugated synthetic peptide corresponding to a region of CHAT gene product shared between the human (P\_066264, NCBI) and mouse (Q03059, NCBI) sequences.

**Host:** Chicken

**Reactivity:** Human, Mouse

**Applications:** ICC, IHC  
(See our web site product page for detailed applications information)

**Protocols:** See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Form:** Liquid

**Purification:** Antigen affinity purification

**Isotype:** IgY

#### Recommend

**Usage:** Immunocytochemistry(1:1000-1:2000)  
Immunohistochemistry(1:1000-1:2000)  
The optimal working dilution should be determined by the end user.

**Storage Buffer:** In PBS, pH 7.2 (0.02% sodium azide)

**Storage Instruction:** Store at 4°C and avoid from light.  
Aliquot to avoid repeated freezing and thawing.

**Entrez GeneID:** 1103|12647

**Gene Symbol:** CHAT

**Gene Alias:** CMS1A, CMS1A2