

Datasheet

CORO1A polyclonal antibody

Catalog Number: PAB29067

Regulatory Status: For research use only (RUO)

Product Description: Chicken polyclonal antibody raised against recombinant Human CORO1A.

Immunogen: A KLH-conjugated synthetic peptide corresponding to a region of CORO1A gene product shared between the mouse (NM_009898) and human (NM_007074) gene products.

Host: Chicken

Reactivity: Human, Mouse

Applications: ICC, IHC
(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Purification: Antigen affinity purification

Isotype: IgY

Recommend

Usage: Immunocytochemistry(1:1000-1:2000)
Immunohistochemistry(1:1000-1:2000)
The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH 7.2 (0.02% sodium azide)

Storage Instruction: Store at 4°C and avoid from light.
Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 11151|12721

Gene Symbol: CORO1A

Gene Alias: CLABP, CLIPINA, FLJ41407, HCORO1, MGC117380, TACO, p57