

Datasheet

DCX polyclonal antibody

Catalog Number: PAB29068

Regulatory Status: For research use only (RUO)

Product Description: Chicken polyclonal antibody raised against recombinant Human DCX.

Immunogen: Two different KLH-conjugated synthetic peptides corresponding to different regions of DCX gene product, shared between the human (CAA06617.1, NCBI) and mouse (AAT58219.1, NCBI) sequences.

Host: Chicken

Reactivity: Human, Mouse

Applications: ICC, IHC
(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Purification: Antigen affinity purification

Isotype: IgY

Recommend

Usage: Immunocytochemistry(1:1000-1:2000)
Immunohistochemistry(1:1000-1:2000)
The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH 7.2 (0.02% sodium azide)

Storage Instruction: Store at 4°C and avoid from light.
Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 1641|13193

Gene Symbol: DCX

Gene Alias: DBCN, DC, LISX, SCLH, XLIS