

Datasheet

GRM2 polyclonal antibody

Catalog Number: PAB29070

Regulatory Status: For research use only (RUO)

Product Description: Chicken polyclonal antibody raised against recombinant Human GRM2.

Immunogen: A KLH-conjugated synthetic peptide corresponding to a unique sequence within GRM2 gene product, shared between the human (Q14416.2, NCBI) and mouse (Q14B12.2, NCBI) sequences.

Host: Chicken

Reactivity: Human, Mouse

Applications: ICC, IHC
(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Purification: Antigen affinity purification

Isotype: IgY

Recommend

Usage: Immunocytochemistry(1:1000-1:2000)

Immunohistochemistry(1:1000-1:2000)

The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH 7.2 (0.02% sodium azide)

Storage Instruction: Store at 4°C and avoid from light. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 2912|108068

Gene Symbol: GRM2

Gene Alias: GLUR2, GPRC1B, MGLUR2, mGlu2