

Datasheet

GRM3 polyclonal antibody

Catalog Number: PAB29071

Regulatory Status: For research use only (RUO)

Product Description: Chicken polyclonal antibody raised against recombinant Human GRM3.

Immunogen: A KLH-conjugated synthetic peptide corresponding to a unique sequence within GRM3 gene product, shared between the human (Q14832.2, NCBI) and mouse (Q9QYS2.1, NCBI) sequences.

Host: Chicken

Reactivity: Human, Mouse

Applications: ICC, IHC
(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Purification: Antigen affinity purification

Isotype: IgY

Recommend

Usage: Immunocytochemistry(1:1000-1:2000)
Immunohistochemistry(1:1000-1:2000)
The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH 7.2 (0.02% sodium azide)

Storage Instruction: Store at 4°C and avoid from light.
Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 2913|108069

Gene Symbol: GRM3

Gene Alias: GLUR3, GPRC1C, MGLUR3, mGlu3