

Datasheet

GRM5 polyclonal antibody

Catalog Number: PAB29072

Regulatory Status: For research use only (RUO)

Product Description: Chicken polyclonal antibody raised against recombinant Human GRM5.

Immunogen: Three different KLH-conjugated synthetic peptides corresponding to different regions of GRM5 gene product, shared between the human (NP_001137303.1, NCBI) and mouse (Q3UVX5.2, NCBI) sequences.

Host: Chicken

Reactivity: Human, Mouse

Applications: ICC, IHC

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Purification: Antigen affinity purification

Isotype: IgY

Recommend

Usage: Immunocytochemistry(1:1000-1:2000)

Immunohistochemistry(1:1000-1:2000)

The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH 7.2 (0.02% sodium azide)

Storage Instruction: Store at 4°C and avoid from light. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 2915|108071

Gene Symbol: GRM5

Gene Alias: GPRC1E, MGLUR5, mGlu5