

Datasheet

GAD1 polyclonal antibody

Catalog Number: PAB29073

Regulatory Status: For research use only (RUO)

Product Description: Chicken polyclonal antibody raised against recombinant Human GAD1.

Immunogen: A KLH-conjugated synthetic peptide corresponding to a region near the C-terminus of GAD1 gene product, and was 100% conserved between the human (Q99259, NCBI), mouse (P48318, NCBI) and rat (NP_058703, NCBI) gene products.

Host: Chicken

Reactivity: Human, Mouse, Rat

Applications: ICC, IHC

(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Purification: Antigen affinity purification

Isotype: IgY

Recommend

Usage: Immunocytochemistry(1:1000-1:2000)

Immunohistochemistry(1:1000-1:2000)

The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH 7.2 (0.02% sodium azide)

Storage Instruction: Store at 4°C and avoid from light. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 2571|14415|24379

Gene Symbol: GAD1

Gene Alias: FLJ45882, GAD, SCP