

## Datasheet

### Gfap polyclonal antibody

**Catalog Number:** PAB29075

**Regulatory Status:** For research use only (RUO)

**Product Description:** Chicken polyclonal antibody raised against recombinant Rat Gfap.

**Immunogen:** *E. coli*-expressed recombinant protein.

**Host:** Chicken

**Reactivity:** Rat

**Applications:** ICC, IHC, WB  
(See our web site product page for detailed applications information)

**Protocols:** See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Form:** Liquid

**Purification:** Antigen affinity purification

**Isotype:** IgY

#### Recommend

**Usage:** Immunocytochemistry(1:2000-1:5000)  
Immunohistochemistry(1:2000-1:5000)  
Western Blot(1:5000-1:10000)  
The optimal working dilution should be determined by the end user.

**Storage Buffer:** In PBS, pH 7.2 (0.02% sodium azide)

**Storage Instruction:** Store at 4° C and avoid from light.  
Aliquot to avoid repeated freezing and thawing.

**Entrez GeneID:** 24387

**Gene Symbol:** Gfap

**Gene Alias:** -