

## Datasheet

### Ilgam polyclonal antibody

**Catalog Number:** PAB29076

**Regulatory Status:** For research use only (RUO)

**Product Description:** Chicken polyclonal antibody raised against recombinant Mouse Ilgam.

**Immunogen:** A KLH-conjugated synthetic peptide corresponding to a region near the center of the extracellular domain of Ilgam gene product.

**Host:** Chicken

**Reactivity:** Mouse

**Applications:** ICC, IHC  
(See our web site product page for detailed applications information)

**Protocols:** See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Form:** Liquid

**Purification:** Antigen affinity purification

**Isotype:** IgY

**Recommend Usage:** Immunocytochemistry  
Immunohistochemistry  
The optimal working dilution should be determined by the end user.

**Storage Buffer:** In PBS, pH 7.2 (0.02% sodium azide)

**Storage Instruction:** Store at 4° C and avoid from light.  
Aliquot to avoid repeated freezing and thawing.

**Entrez GeneID:** 16409

**Gene Symbol:** Ilgam

**Gene Alias:** CD11b/CD18, CR3, CR3A, Cd11b, F730045J24Rik, Ly-40, MAC1, Mac-1, Mac-1a