

Datasheet

MAP2 polyclonal antibody

Catalog Number: PAB29077

Regulatory Status: For research use only (RUO)

Product Description: Chicken polyclonal antibody raised against recombinant Human MAP2.

Immunogen: Two different KLH-conjugated synthetic peptides corresponding to different regions of MAP2 gene product, shared between the human (NP_002365, NCBI) and mouse (P20357, NCBI) sequences.

Host: Chicken

Reactivity: Human, Mouse

Applications: ICC, IHC, WB
(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Purification: Antigen affinity purification

Isotype: IgY

Recommend

Usage: Immunocytochemistry(1:1000-1:2000)

Immunohistochemistry(1:1000-1:2000)

Western Blot(1:2000-1:5000)

The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH 7.2 (0.02% sodium azide)

Storage Instruction: Store at 4°C and avoid from light. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 4133|17756

Gene Symbol: MAP2

Gene Alias: DKFZp686I2148, MAP2A, MAP2B, MAP2C