

Datasheet

MBP polyclonal antibody

Catalog Number: PAB29078

Regulatory Status: For research use only (RUO)

Product Description: Chicken polyclonal antibody raised against recombinant Human MBP.

Immunogen: A KLH-conjugated synthetic peptide shared between the human (NP_002376, NCBI) and mouse (NP_034907, NCBI) sequences.

Host: Chicken

Reactivity: Human, Mouse

Applications: ICC, IHC, WB
(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Purification: Antigen affinity purification

Isotype: IgY

Recommend

Usage: Immunocytochemistry(1:1000-1:2000)

Immunohistochemistry(1:1000-1:2000)

Western Blot(1:2000-1:5000)

The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH 7.2 (0.02% sodium azide)

Storage Instruction: Store at 4°C and avoid from light. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 4155|17196

Gene Symbol: MBP

Gene Alias: MGC99675