

## Datasheet

### Nes polyclonal antibody

**Catalog Number:** PAB29079

**Regulatory Status:** For research use only (RUO)

**Product Description:** Chicken polyclonal antibody raised against recombinant Mouse Nes.

**Immunogen:** Three different KLH-conjugated synthetic peptides corresponding to different regions of the mouse Nes gene product (NP\_057910, NCBI), shared with the rat (NM\_012987, NCBI) protein sequence.

**Host:** Chicken

**Reactivity:** Mouse,Rat

**Applications:** ICC, IHC, WB  
(See our web site product page for detailed applications information)

**Protocols:** See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Form:** Liquid

**Purification:** Antigen affinity purification

**Isotype:** IgY

#### Recommend

**Usage:** Immunocytochemistry(1:10000-1:20000)

Immunohistochemistry(1:10000-1:20000)

Western Blot(1:20000-1:50000)

The optimal working dilution should be determined by the end user.

**Storage Buffer:** In PBS, pH 7.2 (0.02% sodium azide)

**Storage Instruction:** Store at 4°C and avoid from light. Aliquot to avoid repeated freezing and thawing.

**Entrez GeneID:** 18008|25491

**Gene Symbol:** Nes

**Gene Alias:** AA166324, C78523, ESTM46