

Datasheet

Nefm polyclonal antibody

Catalog Number: PAB29083

Regulatory Status: For research use only (RUO)

Product Description: Chicken polyclonal antibody raised against recombinant Mouse Nefm.

Immunogen: Purified Nefm protein (160 kDa) isolated from adult mouse brains.

Host: Chicken

Reactivity: Mouse

Applications: ICC, IHC, WB
(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Purification: IgY Fraction

Isotype: IgY

Recommend

Usage: Immunocytochemistry(1:2000-1:5000)

Immunohistochemistry(1:2000-1:5000)

Western Blot(1:5000-1:10000)

The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH 7.2 (0.2% BSA, 0.02% sodium azide)

Storage Instruction: Store at 4°C and avoid from light. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 18040

Gene Symbol: Nefm

Gene Alias: NF-M, NF160, NF165, Nef3, Nfm