

Datasheet

ENO2 polyclonal antibody

Catalog Number: PAB29084

Regulatory Status: For research use only (RUO)

Product Description: Chicken polyclonal antibody raised against recombinant Human ENO2.

Immunogen: Two different KLH-conjugated synthetic peptides corresponding to different regions of the ENO2 gene product, shared between the human (NP_001966, NCBI) and rat (AAA41119, NCBI) sequences.

Host: Chicken

Reactivity: Human, Rat

Applications: ICC, IHC, WB
(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Purification: Antigen affinity purification

Isotype: IgY

Recommend

Usage: Immunocytochemistry(1:2000-1:5000)

Immunohistochemistry(1:2000-1:5000)

Western Blot(1:5000-1:10000)

The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH 7.2 (0.02% sodium azide)

Storage Instruction: Store at 4°C and avoid from light. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 2026|24334

Gene Symbol: ENO2

Gene Alias: NSE