

Datasheet

hCG_1776007 polyclonal antibody

Gene Alias: FLJ56884, FLJ58356, FOX3, HRNBP3

Catalog Number: PAB29085

Regulatory Status: For research use only (RUO)

Product Description: Chicken polyclonal antibody raised against recombinant Human hCG_1776007.

Immunogen: A KLH-conjugated synthetic peptide corresponding to the hCG_1776007 gene product, shared between the human (NP_001076044, NCBI) and mouse (NP_001020102, NCBI) sequences.

Host: Chicken

Reactivity: Human, Mouse

Applications: ICC, IHC
(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Purification: Antigen affinity purification

Isotype: IgY

Recommend

Usage: Immunocytochemistry(1:1000-1:2000)
Immunohistochemistry(1:1000-1:2000)
The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH 7.2 (50% glycerol, 0.02% sodium azide)

Storage Instruction: Store at 4°C and avoid from light.
Do not freeze unless for long term storage, store at -20°C with 50% glycerol.
Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 146713|52897

Gene Symbol: hCG_1776007