

Datasheet

Acpp polyclonal antibody

Catalog Number: PAB29086

Regulatory Status: For research use only (RUO)

Product Description: Chicken polyclonal antibody raised against recombinant Mouse Acpp.

Immunogen: Recombinant mouse Acpp protein expressed in baculoviral-delivery system.

Host: Chicken

Reactivity: Mouse

Applications: ICC, IHC
(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Purification: IgY Fraction

Isotype: IgY

Recommend Usage: Immunocytochemistry
Immunohistochemistry(1:500-1:1000)
The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH 7.2 (0.02% sodium azide)

Storage Instruction: Store at 4°C and avoid from light.
For long term storage, dilute the antibody with sterile glycerol and store it at -20°C.
Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 56318

Gene Symbol: Acpp

Gene Alias: A030005E02Rik, AI324033, Lap, PAP, Ppal