

## Datasheet

### Prph polyclonal antibody

**Catalog Number:** PAB29087

**Regulatory Status:** For research use only (RUO)

**Product Description:** Chicken polyclonal antibody raised against recombinant Rat Prph.

**Immunogen:** Recombinant full-length rat Prph protein (57 kDa).

**Host:** Chicken

**Reactivity:** Rat

**Applications:** ICC, IHC, WB  
(See our web site product page for detailed applications information)

**Protocols:** See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Form:** Liquid

**Purification:** IgY Fraction

**Isotype:** IgY

#### Recommend

**Usage:** Immunocytochemistry(1:2000-1:5000)

Immunohistochemistry(1:2000-1:5000)

Western Blot(1:5000-1:10000)

The optimal working dilution should be determined by the end user.

**Storage Buffer:** In PBS, pH 7.2 (0.02% sodium azide)

**Storage Instruction:** Store at 4°C and avoid from light.

For long term storage, dilute the antibody with sterile glycerol and store it at -20°C.

Aliquot to avoid repeated freezing and thawing.

**Entrez GeneID:** 24688

**Gene Symbol:** Prph

**Gene Alias:** Perf, Prph1