

Datasheet

PLP1 polyclonal antibody

Catalog Number: PAB29088

Regulatory Status: For research use only (RUO)

Product Description: Chicken polyclonal antibody raised against recombinant Human PLP1.

Immunogen: A KLH-conjugated synthetic peptide corresponding to human and mouse PLP1.(NP_000524, NP_955772, P60202)

Host: Chicken

Reactivity: Human, Mouse

Applications: ICC, IHC
(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Purification: Antigen affinity purification

Isotype: IgY

Recommend

Usage: Immunocytochemistry(1:1000-1:2000)
Immunohistochemistry(1:1000-1:2000)
The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH 7.2 (0.02% sodium azide)

Storage Instruction: Store at 4°C and avoid from light.
For long term storage, dilute the antibody with sterile glycerol and store it at -20°C.
Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 5354|18823

Gene Symbol: PLP1

Gene Alias: HLD1, MMPL, PLP, PLP/DM20, PMD, SPG2