

Datasheet

SYT1 polyclonal antibody

Catalog Number: PAB29091

Regulatory Status: For research use only (RUO)

Product Description: Chicken polyclonal antibody raised against recombinant Human SYT1.

Immunogen: Four different KLH-conjugated synthetic peptides corresponding to different regions of the SYT1 gene product, shared between the human (NP_005630, NCBI) and mouse (NP_033332, NCBI) sequences.

Host: Chicken

Reactivity: Human, Mouse

Applications: ICC, IHC, WB
(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Purification: Antigen affinity purification

Isotype: IgY

Recommend

Usage: Immunocytochemistry(1:1000-1:2000)

Immunohistochemistry(1:1000-1:2000)

Western Blot(1:2000-1:5000)

The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH 7.2 (0.02% sodium azide)

Storage Instruction: Store at 4°C and avoid from light. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 6857|20979

Gene Symbol: SYT1

Gene Alias: DKFZp781D2042, P65, SVP65, SYT