

## Datasheet

### MAPT polyclonal antibody

MGC138549, MSTD, MTBT1, MTBT2, PPND, TAU

**Catalog Number:** PAB29092**Regulatory Status:** For research use only (RUO)**Product Description:** Chicken polyclonal antibody raised against recombinant Human MAPT.**Immunogen:** Two different KLH-conjugated synthetic peptides corresponding to different regions of the MAPT sequences, shared between the mouse and human proteins, and did not contain any of the serine and threonine residues known to undergo phosphorylation.**Host:** Chicken**Reactivity:** Human, Mouse**Applications:** ICC, IHC, WB

(See our web site product page for detailed applications information)

**Protocols:** See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols**Form:** Liquid**Purification:** Antigen affinity purification**Isotype:** IgY**Recommend****Usage:** Immunocytochemistry(1:1000-1:2000)

Immunohistochemistry(1:1000-1:2000)

Western Blot

The optimal working dilution should be determined by the end user.

**Storage Buffer:** In PBS, pH 7.2 (0.02% sodium azide)**Storage Instruction:** Store at 4°C and avoid from light. Aliquot to avoid repeated freezing and thawing.**Entrez GeneID:** 4137|17762**Gene Symbol:** MAPT**Gene Alias:** DDPAC, FLJ31424, FTDP-17, MAPTL,