

## Datasheet

### TUBB3 polyclonal antibody

**Catalog Number:** PAB29093

**Regulatory Status:** For research use only (RUO)

**Product Description:** Chicken polyclonal antibody raised against recombinant Human TUBB3.

**Immunogen:** Three different KLH-conjugated synthetic peptides corresponding to different regions of the TUBB3 gene product, shared between the human (AAL28094, NCBI) and rat (AAM28438, NCBI) sequences.

**Host:** Chicken

**Reactivity:** Human, Rat

**Applications:** ICC, IHC, WB  
(See our web site product page for detailed applications information)

**Protocols:** See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Form:** Liquid

**Purification:** Antigen affinity purification

**Isotype:** IgY

#### Recommend

**Usage:** Immunocytochemistry(1:1000-1:2000)

Immunohistochemistry(1:1000-1:2000)

Western Blot(1:2000-1:5000)

The optimal working dilution should be determined by the end user.

**Storage Buffer:** In PBS, pH 7.2 (0.02% sodium azide)

**Storage Instruction:** Store at 4°C and avoid from light. Aliquot to avoid repeated freezing and thawing.

**Entrez GeneID:** 10381|246118

**Gene Symbol:** TUBB3

**Gene Alias:** MC1R, TUBB4, beta-4