

Datasheet

Goat anti-Chicken IgY secondary antibody (Alkaline Phosphatase)

Catalog Number: PAB29096

Regulatory Status: For research use only (RUO)

Product Description: Goat polyclonal antibody raised against Chicken IgY. The antibody is conjugated with Alkaline Phosphatase.

Immunogen: Purified chicken IgY (whole protein) emulsified in Freund's adjuvant.

Host: Goat

Reactivity: Chicken

Applications: ELISA, WB
(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Conjugation: Alkaline Phosphatase

Purification: Affinity purification

Recommend Usage: ELISA
Western Blot(1:1000)
The optimal working dilution should be determined by the end user.

Storage Buffer: In 50mM HEPES, pH 7.0, 0.1M NaCl, 1.0mM MgCl₂, 0.1mM ZnCl₂ (0.2% BSA, 0.02% sodium azide).

Storage Instruction: Store at 4°C and avoid from light. For long term storage, dilute the antibody with sterile glycerol and store it at -20°C. Aliquot to avoid repeated freezing and thawing.