

## Datasheet

### Goat anti-Chicken IgY secondary antibody (Biotin)

**Catalog Number:** PAB29097

**Regulatory Status:** For research use only (RUO)

**Product Description:** Goat polyclonal antibody raised against Chicken IgY. The antibody is conjugated with Biotin.

**Immunogen:** Purified chicken IgY (whole protein) emulsified in Freund's adjuvant.

**Host:** Goat

**Reactivity:** Chicken

**Applications:** ELISA, ICC, IHC, WB  
(See our web site product page for detailed applications information)

**Protocols:** See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Form:** Liquid

**Conjugation:** Biotin

**Purification:** Affinity purification

**Recommend Usage:** ELISA

Immunocytochemistry

Immunohistochemistry

Western Blot(1:5000)

The optimal working dilution should be determined by the end user.

**Storage Buffer:** In PBS, pH 7.2 (0.02% sodium azide)

**Storage Instruction:** Store at 4°C and avoid from light.  
For long term storage, dilute the antibody with sterile glycerol and store it at -20°C.  
Aliquot to avoid repeated freezing and thawing.