

Datasheet

Chicken anti-Goat IgG secondary antibody (HRP)

Catalog Number: PAB29101

Regulatory Status: For research use only (RUO)

Product Description: Chicken polyclonal antibody raised against Goat IgG. The antibody is conjugated with HRP.

Immunogen: Purified goat IgG (whole protein) emulsified in Freund's adjuvant.

Host: Chicken

Reactivity: Goat

Applications: ELISA, WB
(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Conjugation: HRP

Purification: Antigen affinity purification

Isotype: IgY

Recommend Usage: ELISA
Western Blot(1:5000)
The optimal working dilution should be determined by the end user.

Storage Buffer: In 50mM HEPES, pH 7.0, 0.1M NaCl, (0.2% BSA, 0.01% thimerosal).

Storage Instruction: Store at 4°C and avoid from light.
For long term storage, dilute the antibody with sterile glycerol and store it at -20°C.
Aliquot to avoid repeated freezing and thawing.