

Datasheet

Chicken anti-Human IgG secondary antibody (FITC)

Catalog Number: PAB29103

Regulatory Status: For research use only (RUO)

Product Description: Chicken polyclonal antibody raised against Human IgG. The antibody is conjugated with FITC.

Immunogen: Purified human IgG (whole protein).

Host: Chicken

Reactivity: Human

Applications: ICC, IHC

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Conjugation: FITC

Purification: Antigen affinity purification

Isotype: IgY

Recommend Usage: Immunocytochemistry

Immunohistochemistry(1:250)

The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH 7.2 (0.02% sodium azide)

Storage Instruction: Store at 4°C and avoid from light.

For long term storage, dilute the antibody with sterile glycerol and store it at -20°C.

Aliquot to avoid repeated freezing and thawing.