

Datasheet

Chicken anti-Mouse IgG secondary antibody

Catalog Number: PAB29108

Regulatory Status: For research use only (RUO)

Product Description: Chicken polyclonal antibody raised against Mouse IgG.

Immunogen: Purified mouse IgG (whole protein).

Host: Chicken

Reactivity: Mouse

Applications: ELISA, ICC, IHC, WB
(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Purification: Antigen affinity purification

Isotype: IgY

Recommend Usage: ELISA

Immunocytochemistry

Immunohistochemistry(1:250)

Western Blot(1:1000)

The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH 7.2 (0.02% sodium azide)

Storage Instruction: Store at 4°C and avoid from light.

For long term storage, dilute the antibody with sterile glycerol and store it at -20°C.

Aliquot to avoid repeated freezing and thawing.