

## Datasheet

### Chicken anti-Rabbit IgG secondary antibody (FITC)

**Catalog Number:** PAB29109

**Regulatory Status:** For research use only (RUO)

**Product Description:** Chicken polyclonal antibody raised against Rabbit IgG. The antibody is conjugated with FITC.

**Immunogen:** Purified rabbit IgG (whole protein).

**Host:** Chicken

**Reactivity:** Rabbit

**Applications:** ICC, IHC

(See our web site product page for detailed applications information)

**Protocols:** See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Form:** Liquid

**Conjugation:** FITC

**Purification:** Antigen affinity purification

**Isotype:** IgY

**Recommend Usage:** Immunocytochemistry

Immunohistochemistry(1:250)

The optimal working dilution should be determined by the end user.

**Storage Buffer:** In PBS, pH 7.2 (0.2% BSA, 0.02% sodium azide)

**Storage Instruction:** Store at 4°C and avoid from light.

For long term storage, dilute the antibody with sterile glycerol and store it at -20°C.

Aliquot to avoid repeated freezing and thawing.