

Datasheet

Goat anti-Chicken IgY secondary antibody

Catalog Number: PAB29112

Regulatory Status: For research use only (RUO)

Product Description: Goat polyclonal antibody raised against Chicken IgY.

Immunogen: Purified chicken IgY (whole protein) emulsified in Freund's adjuvant.

Host: Goat

Reactivity: Chicken

Applications: ELISA, ICC, IHC, IP, WB
(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Purification: Affinity purification

Recommend Usage: ELISA
Immunocytochemistry
Immunohistochemistry(1:250)
Immunoprecipitation(1:250)
Western Blot(1:250)

The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH 7.2 (0.02% sodium azide)

Storage Instruction: Store at 4°C and avoid from light.
For long term storage, dilute the antibody with sterile glycerol and store it at -20°C.
Aliquot to avoid repeated freezing and thawing.