

Datasheet

SPP1 polyclonal antibody (FITC)

Catalog Number: PAB29193

Regulatory Status: For research use only (RUO)

Product Description: Chicken polyclonal antibody raised against the synthetic peptide of Human Osteopontin-c with FITC conjugation.

Immunogen: A synthetic peptide corresponding to the splice junction of osteopontin-c.

Sequence: SEEKQNAVSC

Host: Chicken

Reactivity: Human

Applications: IHC
(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Conjugation: FITC

Purification: Antigen affinity purification

Isotype: IgY

Recommend Usage: Immunohistochemistry
The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH 7.2 (0.05% sodium azide, 50% glycerol)

Storage Instruction: Store at 4°C.
Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 6696

Gene Symbol: SPP1

Gene Alias: BNSP, BSPI, ETA-1, MGC110940, OPN