

Datasheet

c-Myc tag polyclonal antibody

Catalog Number: PAB29202

Regulatory Status: For research use only (RUO)

Product Description: Chicken polyclonal antibody raised against c-Myc tag.

Immunogen: Asynthetic peptide corresponding to c-Myc tag and coupled to KLH.

Sequence: EQKLISEEDL

Host: Chicken

Reactivity: Human

Applications: ELISA, Flow Cyt, WB
(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Purification: Antigen affinity purification

Isotype: IgY

Recommend Usage: ELISA (Capture: 1:50-1:200,

Detection: 1:5000-1:20000)

Flow Cytometry (1:100)

Western Blot (1:5000)

The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH 7.2 (0.09% sodium azide)

Storage Instruction: Store at 4°C.

Aliquot to avoid repeated freezing and thawing.