

Datasheet

MKI67 polyclonal antibody

Catalog Number: PAB29307

Regulatory Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against recombinant human MKI67.

Immunogen: Recombinant protein corresponding to amino acids 16-159 of human MKI67.

Sequence:

DGPHFPLSLSTCLFGRGIECDIRIQLPVVSKQHCKIEIH
EQEAILHNFSSTNPTQVNGSVIDEPVRLKHGDVITIIDR
SFRYENESLQSGRKSTEFPRKIREQEPARRVSRSSFS
SDPDEKAQDSKAYSKITEGKVSQNPQVHI

Host: Rabbit

Reactivity: Human

Applications: IF, IHC-P

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Purification: Antigen affinity purification

Isotype: IgG

Recommend Usage: Immunofluorescence (1-4 ug/mL)

Immunohistochemistry (1:50-1:200)

The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH 7.2 (40% glycerol, 0.02% sodium azide)

Storage Instruction: Store at 4°C. For long term storage store at -20°C.

Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 4288

Gene Symbol: MKI67

Gene Alias: KIA, Ki-67

Gene Summary: This gene encodes a nuclear protein that is associated with and may be necessary for cellular proliferation. Alternatively spliced transcript variants have been described. A related pseudogene exists on chromosome X. [provided by RefSeq]