

Datasheet

BrdU polyclonal antibody

Catalog Number: PAB29679

Regulatory Status: For research use only (RUO)

Product Description: Chicken polyclonal antibody raised against BrdU.

Immunogen: BrdU

Host: Chicken

Applications: ELISA, IHC
(See our web site product page for detailed applications information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Purification: Antigen Affinity Purification

Isotype: IgY

Recommend Usage: ELISA
Immunohistochemistry (2-10 ug/mL)
The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH 7.2 (0.09% sodium azide)

Storage Instruction: Store at 4°C.
Aliquot to avoid repeated freezing and thawing.