

Datasheet

Goat Anti-rat IgG Fc secondary antibody

Catalog Number: PAB29754

Regulatory Status: For research use only (RUO)

Product Description: Goat polyclonal antibody raised against rat IgG Fc.

Immunogen: Rat IgG Fc fragment.

Host: Goat

Reactivity: Rat

Applications: ELISA, WB

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: It reacts only with Rat IgG by immunodiffusion and IEP techniques. This does not imply an equal balance of reactivity to all IgG subclasses.

Form: Liquid

Purification: Antigen Affinity Purified

Isotype: IgG

Recommend Usage: ELISA (2 ug/mL)

Western Blot (0.05-0.5 ug/mL)

The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH 7.2 (0.09 % sodium azide)

Storage Instruction: Store at 4° C.