

Datasheet

VERO cells polyclonal antibody (FITC)

Catalog Number: PAB29824

Regulatory Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against VERO cells.

Immunogen: Vero-76 cells (ATCC CRL-1587).

Host: Rabbit

Reactivity: African green monkey

Applications: Flow Cyt, IF, IP, WB
(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Conjugation: FITC

Purification: Protein G purification

Purity: >95%

Isotype: IgG

Recommend Usage: Flow Cytometry

Immunofluorescence

Immunoprecipitation

Western Blot

The optimal working dilution should be determined by the end user.

Storage Buffer: In Phosphate buffered saline, pH 7.4

Storage Instruction: Store at -20°C.
Aliquot to avoid repeated freezing and thawing.