

## Datasheet

### THY1 polyclonal antibody

**Gene Symbol:** THY1

**Catalog Number:** PAB30354

**Gene Alias:** CD90, FLJ33325

**Regulatory Status:** For research use only (RUO)

**Product Description:** Rabbit polyclonal antibody raised against partial recombinant human THY1.

**Immunogen:** Recombinant protein corresponding to human THY1.

**Sequence:**

VTSLTACLVDQSLRLDCRHENTSSSPIQYEFSLTRETK  
KHVLFGTVGVP EHTYRSRTNFTSKYNMKVLYLSAFTS  
KDEGTYTCALHHS GHSPPISSQNVTVLRDKLVKCEGIS  
LLAQNTSWL

**Host:** Rabbit

**Reactivity:** Human

**Applications:** IHC-P

(See our web site product page for detailed applications information)

**Protocols:** See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Form:** Liquid

**Purification:** Antigen affinity purification

**Isotype:** IgG

**Recommend Usage:** Immunohistochemistry

(Formalin/PFA-fixed paraffin-embedded sections)

(1:20-1:50)

The optimal working dilution should be determined by the end user.

**Storage Buffer:** In PBS, pH 7.2 (40% glycerol, 0.02% sodium azide).

**Storage Instruction:** Store at 4°C. For long term storage store at -20°C.

Aliquot to avoid repeated freezing and thawing.

**Entrez GeneID:** 7070