

## Datasheet

### POU5F1 polyclonal antibody

**Catalog Number:** PAB31278

**Regulation Status:** For research use only (RUO)

**Product Description:** Rabbit polyclonal antibody raised against synthetic peptide of POU5F1.

**Immunogen:** A synthetic peptide (conjugated with KLH) corresponding to different parts of human POU5F1.

**Host:** Rabbit

**Reactivity:** Human

**Applications:** ChIP, ChIP-Seq, ELISA, WB-Ce  
(See our web site product page for detailed applications information)

**Protocols:** See our web site at  
<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Form:** Liquid

**Purification:** Affinity purification

**Recommend Usage:** ELISA (1:1000-10000)

Western Blot (1:1000)

ChIP (5 ug/CHIP)

The optimal working dilution should be determined by the end user.

**Storage Buffer:** In PBS (0.05% sodium azide, 0.05% proclin 300).

**Storage Instruction:** Store at -20°C. For long term storage store at -80°C.

Aliquot to avoid repeated freezing and thawing.

**Entrez GeneID:** 5460

**Gene Symbol:** POU5F1

**Gene Alias:** MGC22487, OCT3, OCT4, OTF3, OTF4