

Datasheet

KLF4 polyclonal antibody

Catalog Number: PAB6184

Regulatory Status: For research use only (RUO)

Product Description: Goat polyclonal antibody raised against synthetic peptide of KLF4.

Immunogen: A synthetic peptide corresponding to human KLF4.

Sequence: AVSDALLPSFST-C

Host: Goat

Theoretical MW (kDa): 50.1

Reactivity: Human

Applications: ELISA, WB-Ti

(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Purification: Antigen affinity purification

Concentration: 0.5 mg/mL

Recommend Usage: ELISA (1:64000)

Western Blot (0.1-0.3 ug/mL)

The optimal working dilution should be determined by the end user.

Storage Buffer: In Tris saline, pH 7.3 (0.5% BSA, 0.02% sodium azide)

Storage Instruction: Store at -20°C.

Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 9314

Gene Symbol: KLF4

Gene Alias: EZF, GKLF

References:

1. Identification and characterization of a gene encoding a gut-enriched Kruppel-like factor expressed during growth arrest. Shields JM, Christy RJ, Yang VW. J Biol Chem. 1996 Aug 16;271(33):20009-17.