

## Datasheet

### Rabbit Anti-Chicken IgG (H&L) secondary antibody

**Catalog Number:** PAB9177

**Regulatory Status:** For research use only (RUO)

**Product Description:** Rabbit anti-chicken IgG recognizes chicken IgG whole molecule. This secondary antibody was purified using delipidation and defibrination.

**Immunogen:** Chicken IgG whole molecule

**Host:** Rabbit

**Reactivity:** Chicken

**Applications:** ELISA, IHC, WB  
(See our web site product page for detailed applications information)

**Protocols:** See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Form:** Lyophilized

**Conjugation:** Unconjugated

**Purification:** This product was prepared from monospecific antiserum by a delipidation and defibrination. Assay by immunoelectrophoresis resulted in a single precipitin arc against chicken IgG and chicken serum.

**Recommend Usage:** The optimal working dilution should be determined by the end user.

**Storage Buffer:** Lyophilized from 10 mM  $\text{KH}_2\text{PO}_4$ , 150 mM NaCl, pH 7.2

**Storage Instruction:** Store at 4°C on dry atmosphere. After reconstitution with 2 mL of deionized water, store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.